


## ComfyCuff®

|                      |                |   |                    |         |   |
|----------------------|----------------|---|--------------------|---------|---|
| <b>Media Type</b>    | Wristband Cuff | ✓ | <b>Fastening</b>   | Velcro  | ✓ |
| <b>Material Type</b> | Paper          |   | <b>Environment</b> | Indoor  | ✓ |
|                      | Synthetic      | ✓ |                    | Outdoor |   |

### Additional Features

- Super soft foam, padded cuff for superior comfort when identifying patients with delicate skin
- Velcro® fastening means ComfyCuff can be adjusted for optimum comfort, for example as baby's weight fluctuates or if patient suffers from swelling or oedema
- Designed to be used in conjunction with other Z-Band® wristbands as well as LaserBand® wristbands and ankle ID tags for babies
- Can be used in conjunction with a swaddle band and ankle ID tag, so that baby can remain swaddled whilst having their identity confirmed
- Available in sizes suitable for baby/neonate and adult patients



### Suggested Applications

- Patient identification, when used with other printed wristbands or ankle ID tags

### Technical Specifications

|                  | Description  |
|------------------|--|
| <b>Facestock</b> | Super soft polyester and nylon padded cuff                                 |
| <b>Adhesive</b>  | Nylon with double sided tape section to secure printed wristband or ID tag |

**Recommended Storage Conditions:**  
Storage of product before use

2 years duration when stored at 18°C to 23°C at 45% to 55% RH

**Expected Life Span in Application:**  
Following correct application and appropriate dwell time (usually 24hrs) we expect, but do not warrant, a life span as indicated

Indoor use, short term only

### Product Performance and Suitability

The information contained in this document is to be used for guidance only and is not intended for use in setting specifications. All purchasers of Zebra products shall be solely responsible for independently determining if the product conforms to all requirements of their unique application.

For testing of this product, please contact your Supplies or Healthcare Specialist Representative.