



SEE MORE. DO MORE.

Zebra Validation Program Application Test Overview Mother's Milk



Summary

- [Mother's Milk](#) was launched by [Rapid Consulting Services](#) as a newborn care solution developed for neonatal intensive care units in hospitals.
- Rapid Consulting Services saw the need to integrate mobile barcode printing into their software. The Zebra Global Validation team worked with engineers at Rapid Consulting to test and validate Mother's Milk with Zebra's QLn220 mobile printer.



Product Details

Description: Mother's Milk allows monitoring milk feeding in units of neonatal intensive care in Hospitals. Key features include:

- Scan baby's wristband into the application.
- Print matching labels and apply to infant milk bottles for the mother to take home.
- After the bottles have been filled and returned, scan the labels into the application.
- To complete a feeding, scan the baby's wristband to verify a match between bottle and baby.

Verticals: Healthcare

Validated Application Details

ISV Application:	Mother's Milk
Company:	Rapid Consulting Services
Software Version:	2.2.4
Hardware:	Hand-held client, iPhone 5, iPod Touch 5G, Infinite Peripherals LP5 (2D scanner) for iPod or iPhone
Operating System(s):	iOS 5.0
Command Language:	CPCL / Line Print
Formatting Method:	Zebra Link-OS SDK V2.0.214

Zebra Hardware:	Zebra QLn220
Part #:	QN2-AUCA0E00-00
Description:	QLn220 mobile Bluetooth label printer.
Firmware Version:	V68.19.4z
Website:	www.zebra.com/QLn




Zebra Validation Program

Application Test Overview

Mother's Milk



System Diagram





SEE MORE. DO MORE.

Zebra Validation Program Application Test Overview *Mother's Milk*



Validation Certificate



VALIDATED

Rapid Consulting Services
Mother's Milk, 2.2.4
Validated with Zebra Technologies
QLn220
on December 2013



SEE MORE. DO MORE.



Daniel Quagliana
Global Alliance & ISV Engineering Manager
December 31, 2013

