## **COVID-19 Drive-Through Specimen Collection Station Solutions**



Challenge

Optimize COVID-19 drive-through specimen collection with safety, speed and accuracy. **Solution** 

Leverage mobile barcode technology for data collection, positive patient ID, accurate labeling and printing.

Good

## **Laptop-Based/Electronic Health Record Configuration**

Approx 5-6 Minutes | No Physical Contact | Additional Footsteps | Requires Internet Connection







ID scanned, patient data captured



Healthcare technician leaves patient to enter additional data into laptop



Labels printed



Healthcare technician returns to patient



Swab taken



Patient exits

Plus

**Emergency Wristband Kit** 

**Better** 

## **Turn-Key Mobile Configuration**

Under 3 Minutes | No Physical Contact | No Manual Entry | No Internet Connection



Patients line up



ID scanned, data captured



Labels printed



ID scanned, data captured



Barcode Wristband ID printed



Patient immediately taken to hospital care

Swab taken



Patient exits